

ПАРАМЕТРЫ торических контактных линз ACUVUE®

Сила сферы	от plano до -6,00 D (шаг 0,25 D)					от -6,50 D до -9,00 D (шаг 0,5 D)					от +0,25 D до +4,00 D (шаг 0,25 D)			от +4,25 D до +6,00 D (шаг 0,25 D)		
Сила цилиндра	cyl -0,75	cyl -1,25	cyl -1,75	cyl -2,25	cyl -2,75	cyl -0,75	cyl -1,25	cyl -1,75	cyl -2,25	cyl -2,75	cyl -0,75	cyl -1,25	cyl -1,75	cyl -0,75	cyl -1,25	cyl -1,75

10°																
-----	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--

30°				-	-			-	-	-	-	-	-	-	-	-
40°				-	-			-	-	-	-	-	-	-	-	-
50°				-	-			-	-	-	-	-	-	-	-	-
60°				-	-			-	-	-	-	-	-	-	-	-
70°								-	-	-						
80°								-	-	-						
90°								-	-	-						
100°								-	-	-						
110°								-	-	-						
120°				-	-			-	-	-	-	-	-	-	-	-
130°				-	-			-	-	-	-	-	-	-	-	-
140°				-	-			-	-	-	-	-	-	-	-	-
150°				-	-			-	-	-	-	-	-	-	-	-
160°																
170°																
180°																

ACUVUE® OASYS¹ for ASTIGMATISM² with³ HYDRACLEAR® PLUS⁴

ACUVUE® OASYS¹ 1-DAY⁵ with³ HydraLuxe® for ASTIGMATISM^{2E}

1-DAY ACUVUE® MOIST⁶ for ASTIGMATISM² with³ LACREON®



ОНЛАЙН-калькулятор для расчёта оптической силы торических линз ACUVUE® поможет выбрать параметры контактных линз быстро и легко.

Отсканируйте QR-код

